


DATASHEET

Theora Connect Wearable and Theora Link Smartphone App

Keeping families connected with the stylish cellular wearable for active older adults



Theora Link™ Smartphone App
Caregiver app with safe zones and wandering alerts



Theora Connect™ wearable
with selectable watch faces

The Theora Connect™ wearable for older adults, along with the powerful Theora Link™ smartphone app for families and caregivers, provides a complete solution for communication, remotely monitoring activity, and receiving wandering alerts.

The Theora Connect wearable is a lightweight, stylish wristwatch that people enjoy wearing. It's flexible design offers selectable analog, digital, and date/time watch faces. An easy to locate SOS button directly calls a trusted family member in case of emergency, and family can initiate calls quickly if they receive an alert, or just to chat.

Quik Connect™: Auto-Answering Audio on the Theora Connect Wearable

Theora Connect's touchless audio feature enables family members to communicate with loved ones through a speaker on the Theora Connect wearable. No button pushing is required by the wearer to hear your reassuring voice. At home, in assisted living or memory care – across town or across country – families are staying connected!

Safe Zones: Be Notified of Wandering Before it Escalates

In the event of a wandering incident, the Safe Zone feature of the Theora Link application notifies the primary family member on their smartphone. They can easily track the wearer using the built-in GPS, and call the wearer to help and reassure them help is on the way.

Weather App: Know what's happening outside for better planning

Going for a walk, or just want to know if you need a sweater? Our weather app provides an easy to read and graphical display of current weather conditions. Simply tap on the weather icon on the Theora Connect watch face menu to see current temperature, air quality, weather conditions for your current city location. No need for an additional cell phone to check the weather or make and receive calls.



Voice First Technology

A voice prompt confirms the battery is charging. Know your watch is fully charged and ready to go with audio feedback in addition to visual feedback. This hands-free feature is beneficial to those with physical or cognitive limitations or those who prefer audio.

Theora® customers use our stylish cellular wearable to stay active and independent, connect with family, and monitor general wellbeing. Instantly check the weather, count your steps, and only receive calls from trusted contacts. Even includes an SOS button just in case. Enjoy more peace of mind by staying connected and communicating - it's the only device you need to carry!



- ✔ Enables older loved ones to remain active, independent and in their home of choice longer
- ✔ Automatically connect with loved ones in assisted living or care homes, without care staff assistance
- ✔ Stylish, lightweight wristwatch design
- ✔ Proven solution for Alzheimer's, dementia, Autism, and other conditions
- ✔ More peace of mind for busy family members and caregivers
- ✔ Automatic text notifications if loved one should depart from their usual locations
- ✔ Quickly find and communicate with a loved one who's wandered, before it becomes dangerous
- ✔ Theora Link is a free app for family and trusted contacts to keep everyone informed
- ✔ Built-in pedometer measures daily steps for active adults
- ✔ A programmable lockdown mode for the watch face prevents wearers with cognitive issues from inadvertently changing the settings, potentially causing confusion. SOS button and incoming calls remain enabled
- ✔ Investment protection with latest 4G/LTE cellular
- ✔ Optional accessories include lockable clasp, instant wandering alert and audio door alarm

Theora Connect Product Specifications

Short Wristband	5 ¼" to 7" (Optional extender adds another ½" to either band)
Long Wristband	6 ¼" to 7 ¾"
Watch face	1 ¼" diameter
Water resistant	Yes, IP67

Theora Link Smartphone App Specifications

Theora Link supports the Apple® smartphone iOS and Samsung® Android™ OS. Just search for "TheoraLink" via Google Play™ or Apple® App Store to verify your phone's compatibility.



Need Support?

Email: support@theoracare.com
 Short, simple instructional videos at www.theoracare.com/support
 Call: 512-371-6164

